

TABLE OF CONTENTS

ANTIMICROBIAL RESISTANCE GENE CONTROLS	TYPES	PAGE NUMBERS
SHEET 1	PCR controls for ESBL & pAmpC	2-7
SHEET 2	PCR controls for PMQR genes	8-12
SHEET 3	PCR controls for carbapenemase	13-19

SHEET 1

PCR controls for ESBL & pAmpC

Escherichia coli (Biosafety Level: 2)

Pseudomonas aeruginosa (Biosafety Level: 2)

These multidrug-resistant isolates were chosen to represent different variants of Class A (Extended Spectrum β -lactamase, ESBL) & Class C β -lactamases (Plasmid-mediated AmpC, pAmpC) available at this repository.

Revival of Culture from Stab:

1. Open the stab culture tube under a Laminar Air Flow (LAF) cabinet.
2. Using a sterile needle (straight or loop type - either disposable plastic or sterilized reusable metal), stab into the existing line of growth within the medium.
3. Once growth material adheres to the needle, inoculate it into a general or enriched broth medium (e.g., LB or TSB).
4. Incubate the broth culture at 37°C overnight (O/N).
5. On the next day, using a sterile loop, streak the broth culture onto a non-selective agar plate to obtain isolated colonies.
6. Incubate the plate at 37°C overnight to obtain single colonies.

Preparation of 20% Glycerol Stock (for Long-term Preservation):

1. Prepare 50% glycerol by mixing sterile distilled water and glycerol in a 1:1 ratio in a conical flask.
2. Autoclave the mixture at 121°C for 15 minutes and allow it to cool.
3. Dispense 400 µL of the sterile 50% glycerol into each sterile cryovial.
4. Pick colonies from a fresh nutrient agar plate and inoculate into TSB broth to prepare a dense bacterial suspension.
5. Incubate the broth culture at 37°C in a shaker incubator overnight.
6. After incubation, mix 600 µL of the overnight TSB culture with 400 µL of 50% glycerol to obtain a final glycerol concentration of 20%.
7. Mix thoroughly and store at -80°C for long-term preservation.

LAB ID	Organism	Class A β -lactamases (ESBL)	Susceptibility (MIC)	Accession number
AMRIR00300	<i>Klebsiella pneumoniae</i>	<i>bla</i> _{CTX-M-15} <i>bla</i> _{SHV-182} <i>bla</i> _{TEM-1B}	Ampicillin (≥ 32 mg/L, R) Amoxicillin/Clavulanic acid (≥ 32 mg/L, R) Piperacillin/Tazobactam (≥ 128 mg/L, R) Cefuroxime (≥ 64 mg/L, R) Ceftriaxone (≥ 64 mg/L, R) Cefepime (≥ 64 mg/L, R) Ertapenem (≥ 8 mg/L, R) Imipenem (2 mg/L, I) Meropenem (≥ 16 mg/L, R) Amikacin (≤ 2 mg/L, S) Gentamicin (≥ 16 mg/L, R) Nalidixic acid (≥ 32 mg/L, R) Ciprofloxacin (≥ 4 mg/L, R) Colistin (≤ 0.125 mg/L, I) Trimethoprim/Sulfamethoxazole (≥ 320 mg/L, R)	JAWVAF000000000

AMRIR00496	<i>Klebsiella pneumoniae</i>	<i>bla</i> _{CTX-M-15} <i>bla</i> _{SHV-106, 28} <i>bla</i> _{TEM-1A}	Amoxicillin/Clavulanic acid (≥32 mg/L, R) Piperacillin/Tazobactam (≥128 mg/L, R) Cefuroxime (≥64 mg/L, R) Ceftriaxone (≥64 mg/L, R) Cefepime (≥32 mg/L, R) Ertapenem (2 mg/L, R) Imipenem (≤0.25 mg/L, S) Meropenem (≤0.25 mg/L, S) Amikacin (32 mg/L, I) Gentamicin (≥16 mg/L, R) Ciprofloxacin (≥4 mg/L, R) Fosfomycin (32 mg/L, S) Colistin (2 mg/L, I) Trimethoprim/Sulfamethoxazole (≥320 mg/L, R)	JAWVAK000000000
AMRIR00637	<i>Pseudomonas aeruginosa</i>	<i>bla</i> _{VEB-5}	Piperacillin/Tazobactam (≥128 mg/L, R) Ceftazidime (≥64 mg/L, R) Cefepime (≥32 mg/L, R) Aztreonam (≥64 mg/L, R) Imipenem (≥16 mg/L, R) Meropenem (≥16 mg/L, R) Amikacin (≥64 mg/L, R) Ciprofloxacin (≥4 mg/L, R) Levofloxacin (≥8 mg/L, R) Colistin (1 mg/L, I)	JBLAVS000000000
AMRIR00636	<i>Pseudomonas aeruginosa</i>	<i>bla</i> _{VEB-1}	Piperacillin/Tazobactam (64 mg/L, R) Ceftazidime (≥64 mg/L, R) Cefepime (≥32 mg/L, R)	JBLAVR000000000

			Aztreonam (≥ 64 mg/L, R) Imipenem (≥ 16 mg/L, R) Meropenem (≥ 16 mg/L, R) Amikacin (≥ 64 mg/L, R) Ciprofloxacin (≥ 4 mg/L, R) Levofloxacin (≥ 8 mg/L, R) Colistin (1 mg/L, I)	
LAB ID	Organism	Class C β -lactamases (pAmpC)	Susceptibility (MIC)	Accession number
AMRIR00274	<i>Escherichia coli</i>	<i>bla</i> _{DHA-1}	Ampicillin (≥ 32 mg/L, R) Amoxicillin/Clavulanic acid (16 mg/L, I) Piperacillin/Tazobactam (≥ 128 mg/L, R) Cefuroxime (≥ 64 mg/L, R) Ceftriaxone (≥ 64 mg/L, R) Cefepime (32 mg/L, R) Ertapenem (≤ 0.5 mg/L, S) Imipenem (≤ 0.25 mg/L, S) Meropenem (4 mg/L, S) Amikacin (4 mg/L, S) Gentamicin (≤ 1 mg/L, S) Nalidixic acid (≥ 32 mg/L, R) Ciprofloxacin (> 4 mg/L, R) Nitrofurantoin (≤ 16 mg/L, S) Colistin (0.5 mg/L, I) Trimethoprim/Sulfamethoxazole (≥ 320 mg/L, R)	JBLFCW000000000
AMRIR00482	<i>Escherichia coli</i>	<i>bla</i> _{DHA-1}	Amoxicillin/Clavulanic acid (≥ 32 mg/L, R) Piperacillin/Tazobactam (≤ 4 mg/L, S) Cefuroxime (16 mg/L, I)	JBAGCE000000000

			Ceftriaxone (≤ 0.25 mg/L, S) Cefepime (≤ 0.12 mg/L, S) Ertapenem (≤ 0.12 mg/L, S) Imipenem (≤ 0.25 mg/L, S) Meropenem (≤ 0.25 mg/L, S) Amikacin (4 mg/L, S) Gentamicin (≤ 1 mg/L, S) Ciprofloxacin (0.5 mg/L, I) Fosfomycin (≤ 16 mg/L, S) Colistin (0.25 mg/L, I) Trimethoprim/Sulfamethoxazole (≥ 320 mg/L, R)	
AMRIR00460	<i>Escherichia coli</i>	<i>bla</i> _{CMY-42}	Amoxicillin/Clavulanic acid (≥ 32 mg/L, R) Piperacillin/Tazobactam (≥ 128 mg/L, R) Cefuroxime (≥ 64 mg/L, R) Ceftriaxone (≥ 64 mg/L, R) Cefepime (≥ 32 mg/L, R) Ertapenem (0.25 mg/L, S) Imipenem (0.5 mg/L, S) Meropenem (≤ 0.25 mg/L, S) Amikacin (4 mg/L, S) Gentamicin (≤ 1 mg/L, S) Ciprofloxacin (≥ 4 mg/L, R) Fosfomycin (≤ 16 mg/L, S) Colistin (0.5 mg/L, I) Trimethoprim/Sulfamethoxazole (≥ 320 mg/L, R)	JAZEUD000000000
AMRIR00966	<i>Escherichia coli</i>	<i>bla</i> _{CMY-2}	Amoxicillin/Clavulanic acid (≥ 32 mg/L, R) Piperacillin/Tazobactam (≥ 128 mg/L, R) Cefuroxime (≥ 64 mg/L, R)	JBLFDE000000000

			Ceftriaxone (≥ 64 mg/L, R) Cefepime (≥ 32 mg/L, R) Ertapenem (4 mg/L, R) Imipenem (1 mg/L, S) Meropenem (0.5 mg/L, S) Amikacin (4 mg/L, S) Gentamicin (≤ 1 mg/L, S) Ciprofloxacin (≥ 4 mg/L, R) Fosfomycin (≤ 16 mg/L, S) Colistin (0.5 mg/L, I) Trimethoprim/Sulfamethoxazole (≤ 20 mg/L, S)	
--	--	--	---	--

SHEET 2

PCR controls for PMQR genes

Escherichia coli (Biosafety Level: 2)

Klebsiella pneumoniae (Biosafety Level: 2)

These isolates were chosen to represent different variants of Plasmid-mediated Quinolone Resistance Determinants (PMQR) available at the repository.

Revival of Culture from Stab:

8. Open the stab culture tube under a Laminar Air Flow (LAF) cabinet.
9. Using a sterile needle (straight or loop type - either disposable plastic or sterilized reusable metal), stab into the existing line of growth within the medium.
10. Once growth material adheres to the needle, inoculate it into a general or enriched broth medium (e.g., LB or TSB).
11. Incubate the broth culture at 37°C overnight (O/N).
12. On the next day, using a sterile loop, streak the broth culture onto a non-selective agar plate to obtain isolated colonies.
13. Incubate the plate at 37°C overnight to obtain single colonies.

Preparation of 20% Glycerol Stock (for Long-term Preservation):

14. Prepare 50% glycerol by mixing sterile distilled water and glycerol in a 1:1 ratio in a conical flask.
15. Autoclave the mixture at 121°C for 15 minutes and allow it to cool.
16. Dispense 400 µL of the sterile 50% glycerol into each sterile cryovial.
17. Pick colonies from a fresh nutrient agar plate and inoculate into TSB broth to prepare a dense bacterial suspension.
18. Incubate the broth culture at 37°C in a shaker incubator overnight.
19. After incubation, mix 600 µL of the overnight TSB culture with 400 µL of 50% glycerol to obtain a final glycerol concentration of 20%.
20. Mix thoroughly and store at -80°C for long-term preservation.

LAB ID	Organism	Plasmid-mediated Quinolone Resistance Determinants (PMQR)	Susceptibility (MIC)	Accession numbers
AMRIR00460	<i>Escherichia coli</i>	<i>aac(6')-Ib-cr</i>	Amoxicillin/Clavulanic acid (≥ 32 mg/L, R) Piperacillin/Tazobactam (≥ 128 mg/L, R) Cefuroxime (≥ 64 mg/L, R) Ceftriaxone (≥ 64 mg/L, R) Cefepime (≥ 32 mg/L, R) Ertapenem (0.25 mg/L, S) Imipenem (0.5 mg/L, S) Meropenem (≤ 0.25 mg/L, S) Amikacin (4 mg/L, S) Gentamicin (≤ 1 mg/L, S) Ciprofloxacin (≥ 4 mg/L, R) Fosfomycin (≤ 16 mg/L, S)	JAZEUD000000000

			Colistin (0.5 mg/L, I) Trimethoprim/Sulfamethoxazole (>=320 mg/L, R)	
AMRIR00300	<i>Klebsiella pneumoniae</i>	<i>aac(6')-Ib-cr</i>	Ampicillin (>=32 mg/L, R) Amoxicillin/Clavulanic acid (>=32 mg/L, R) Piperacillin/Tazobactam (>=128 mg/L, R) Cefuroxime (>=64 mg/L, R) Ceftriaxone (>=64 mg/L, R) Cefepime (>=64 mg/L, R) Ertapenem (>=8 mg/L, R) Imipenem (2 mg/L, I) Meropenem (>=16 mg/L, R) Amikacin (<=2 mg/L, S) Gentamicin (>=16 mg/L, R) Nalidixic acid (>=32 mg/L, R) Ciprofloxacin (>=4 mg/L, R) Colistin (<=0.125 mg/L, I) Trimethoprim/Sulfamethoxazole (>=320 mg/L, R)	JAWVAF000000000
AMRIR00300	<i>Klebsiella pneumoniae</i>	<i>qnrB1</i>	Ampicillin (>=32 mg/L, R) Amoxicillin/Clavulanic acid (>=32 mg/L, R) Piperacillin/Tazobactam (>=128 mg/L, R) Cefuroxime (>=64 mg/L, R) Ceftriaxone (>=64 mg/L, R) Cefepime (>=64 mg/L, R) Ertapenem (>=8 mg/L, R) Imipenem (2 mg/L, I) Meropenem (>=16 mg/L, R) Amikacin (<=2 mg/L, S) Gentamicin (>=16 mg/L, R)	JAWVAF000000000

			Nalidixic acid (≥ 32 mg/L, R) Ciprofloxacin (≥ 4 mg/L, R) Colistin (≤ 0.125 mg/L, I) Trimethoprim/Sulfamethoxazole (≥ 320 mg/L, R)	
AMRIR00482	<i>Escherichia coli</i>	<i>qnrB4</i>	Amoxicillin/Clavulanic acid (≥ 32 mg/L, R) Piperacillin/Tazobactam (≤ 4 mg/L, S) Cefuroxime (16 mg/L, I) Ceftriaxone (≤ 0.25 mg/L, S) Cefepime (≤ 0.12 mg/L, S) Ertapenem (≤ 0.12 mg/L, S) Imipenem (≤ 0.25 mg/L, S) Meropenem (≤ 0.25 mg/L, S) Amikacin (4 mg/L, S) Gentamicin (≤ 1 mg/L, S) Ciprofloxacin (0.5 mg/L, I) Fosfomycin (≤ 16 mg/L, S) Colistin (0.25 mg/L, I) Trimethoprim/Sulfamethoxazole (≥ 320 mg/L, R)	JBAGCE000000000
AMRIR00300	<i>Klebsiella pneumoniae</i>	<i>qnrS1</i>	Ampicillin (≥ 32 mg/L, R) Amoxicillin/Clavulanic acid (≥ 32 mg/L, R) Piperacillin/Tazobactam (≥ 128 mg/L, R) Cefuroxime (≥ 64 mg/L, R) Ceftriaxone (≥ 64 mg/L, R) Cefepime (≥ 64 mg/L, R) Ertapenem (≥ 8 mg/L, R) Imipenem (2 mg/L, I) Meropenem (≥ 16 mg/L, R) Amikacin (≤ 2 mg/L, S)	JAWVAF000000000

			Gentamicin (≥ 16 mg/L, R) Nalidixic acid (≥ 32 mg/L, R) Ciprofloxacin (≥ 4 mg/L, R) Colistin (≤ 0.125 mg/L, I) Trimethoprim/Sulfamethoxazole (≥ 320 mg/L, R)	
AMRIR00499	<i>Klebsiella pneumoniae</i>	<i>qnrS1</i>	Amoxicillin/Clavulanic acid (≥ 32 mg/L, R) Piperacillin/Tazobactam (≥ 128 mg/L, R) Cefuroxime (≥ 64 mg/L, R) Ceftriaxone (≥ 64 mg/L, R) Cefepime (≥ 32 mg/L, R) Ertapenem (≥ 8 mg/L, R) Imipenem (4 mg/L, R) Meropenem (≥ 16 mg/L, R) Amikacin (32 mg/L, I) Gentamicin (≥ 16 mg/L, R) Ciprofloxacin (≥ 4 mg/L, R) Fosfomycin (32 mg/L, S) Colistin (2 mg/L, I) Trimethoprim/Sulfamethoxazole (≥ 320 mg/L, R)	JAWVAN00000000

SHEET 3

PCR controls for carbapenemase

Escherichia coli (Biosafety Level: 2)

Klebsiella pneumoniae (Biosafety Level: 2)

Pseudomonas aeruginosa (Biosafety Level: 2)

Acinetobacter baumannii (Biosafety Level: 2)

These multidrug-resistant isolates were chosen to represent different variants of carbapenemases available at this repository.

Revival of Culture from Stab:

21. Open the stab culture tube under a Laminar Air Flow (LAF) cabinet.
22. Using a sterile needle (straight or loop type - either disposable plastic or sterilized reusable metal), stab into the existing line of growth within the medium.
23. Once growth material adheres to the needle, inoculate it into a general or enriched broth medium (e.g., LB or TSB).
24. Incubate the broth culture at 37°C overnight (O/N).
25. On the next day, using a sterile loop, streak the broth culture onto a non-selective agar plate to obtain isolated colonies.
26. Incubate the plate at 37°C overnight to obtain single colonies.

Preparation of 20% Glycerol Stock (for Long-term Preservation):

27. Prepare 50% glycerol by mixing sterile distilled water and glycerol in a 1:1 ratio in a conical flask.
28. Autoclave the mixture at 121°C for 15 minutes and allow it to cool.
29. Dispense 400 µL of the sterile 50% glycerol into each sterile cryovial.
30. Pick colonies from a fresh nutrient agar plate and inoculate into TSB broth to prepare a dense bacterial suspension.
31. Incubate the broth culture at 37°C in a shaker incubator overnight.
32. After incubation, mix 600 µL of the overnight TSB culture with 400 µL of 50% glycerol to obtain a final glycerol concentration of 20%.
33. Mix thoroughly and store at -80°C for long-term preservation.

LAB ID	Organism	Class B carbapenemase (Metallo-β-lactamase, MBL)	Susceptibility (MIC)	Accession number
AMRIR00274	<i>Escherichia coli</i>	<i>bla</i> _{NDM-1}	Ampicillin (≥32 mg/L, R) Amoxicillin/Clavulanic acid (16 mg/L, I) Piperacillin/Tazobactam (≥128 mg/L, R) Cefuroxime (≥64 mg/L, R) Ceftriaxone (≥64 mg/L, R) Cefepime (32 mg/L, R) Ertapenem (≤0.5 mg/L, S) Imipenem (≤0.25 mg/L, S) Meropenem (4 mg/L, S) Amikacin (4 mg/L, S) Gentamicin (≤1 mg/L, S) Nalidixic acid (≥32 mg/L, R) Ciprofloxacin (>4 mg/L, R)	JBLFCW000000000

			Nitrofurantoin (<=16 mg/L, S) Colistin (0.5 mg/L, I) Trimethoprim/Sulfamethoxazole (>=320 mg/L, R)	
AMRIR00139	<i>Pseudomonas aeruginosa</i>	<i>bla</i> _{NDM-5}	Ticarcillin/Clavulanic acid (>=128 mg/L, R) Ceftazidime (>=64 mg/L, R) Cefepime (>=64 mg/L, R) Aztreonam (>=64 mg/L, R) Doripenem (>=8 mg/L, R) Imipenem (>=16 mg/L, R) Meropenem (>=16 mg/L, R) Amikacin (>=64 mg/L, R) Gentamicin (>=16 mg/L, R) Ciprofloxacin (>=4 mg/L, R) Levofloxacin (>=8 mg/L, R) Colistin (2 mg/L, I)	JBBBSP000000000
AMRIR00139	<i>Pseudomonas aeruginosa</i>	<i>bla</i> _{VIM-5}	Ticarcillin/Clavulanic acid (>=128 mg/L, R) Ceftazidime (>=64 mg/L, R) Cefepime (>=64 mg/L, R) Aztreonam (>=64 mg/L, R) Doripenem (>=8 mg/L, R) Imipenem (>=16 mg/L, R) Meropenem (>=16 mg/L, R) Amikacin (>=64 mg/L, R) Gentamicin (>=16 mg/L, R) Ciprofloxacin (>=4 mg/L, R) Levofloxacin (>=8 mg/L, R) Colistin (2 mg/L, I)	JBBBSP000000000
AMRIR00154	<i>Pseudomonas</i>	<i>bla</i> _{VIM-2}	Ticarcillin/Clavulanic acid (>=128 mg/L, R)	JBBBJZ000000000

	<i>aeruginosa</i>		Ceftazidime (≥ 64 mg/L, R) Cefepime (≥ 64 mg/L, R) Aztreonam (16 mg/L, I) Doripenem (≥ 8 mg/L, R) Imipenem (≥ 16 mg/L, R) Meropenem (≥ 16 mg/L, R) Amikacin (≥ 64 mg/L, R) Gentamicin (≥ 16 mg/L, R) Ciprofloxacin (≥ 4 mg/L, R) Levofloxacin (≥ 8 mg/L, R) Colistin (1 mg/L, I)	
AMRIR00144	<i>Pseudomonas aeruginosa</i>	<i>bla_{IMP-1}</i>	Ticarcillin/Clavulanic acid (≥ 128 mg/L, R) Piperacillin/Tazobactam (32mg/L, I) Ceftazidime (≥ 64 mg/L, R) Cefepime (≥ 64 mg/L, R) Aztreonam (≥ 64 mg/L, R) Doripenem (≥ 8 mg/L, R) Imipenem (≥ 16 mg/L, R) Meropenem (≥ 16 mg/L, R) Amikacin (32 mg/L, I) Gentamicin (≥ 16 mg/L, R) Ciprofloxacin (≥ 4 mg/L, R) Levofloxacin (≥ 8 mg/L, R) Colistin (1 mg/L, I)	Sanger sequencing
AMRIR00159	<i>Pseudomonas aeruginosa</i>	<i>bla_{IMP-1}</i>	Ticarcillin/Clavulanic acid (≥ 128 mg/L, R) Piperacillin/Tazobactam (64mg/L, I) Ceftazidime (≥ 64 mg/L, R) Cefepime (≥ 64 mg/L, R) Aztreonam (≥ 64 mg/L, R) Doripenem (≥ 8 mg/L, R)	Sanger sequencing

			Imipenem (≥ 16 mg/L, R) Meropenem (≥ 16 mg/L, R) Amikacin (32 mg/L, I) Gentamicin (≥ 16 mg/L, R) Ciprofloxacin (≥ 4 mg/L, R) Levofloxacin (≥ 8 mg/L, R) Colistin (1 mg/L, I)	
AMRIR00111	<i>Pseudomonas aeruginosa</i>	<i>bla</i> _{DIM-1}	Ticarcillin/Clavulanic acid (≥ 128 mg/L, R) Ceftazidime (≥ 64 mg/L, R) Cefepime (≥ 64 mg/L, R) Aztreonam (≥ 64 mg/L, R) Doripenem (≥ 8 mg/L, R) Imipenem (≥ 16 mg/L, R) Meropenem (≥ 16 mg/L, R) Amikacin (≥ 64 mg/L, R) Gentamicin (≥ 16 mg/L, R) Ciprofloxacin (≥ 4 mg/L, R) Levofloxacin (≥ 8 mg/L, R) Colistin (2 mg/L, I)	JBBBSN000000000
LAB ID	Organism	Class A carbapenemase	Susceptibility (MIC)	Accession number
AMRIR00485	<i>Klebsiella pneumoniae</i>	<i>bla</i> _{KPC-2}	Amoxicillin/Clavulanic acid (≥ 32 mg/L, R) Piperacillin/Tazobactam (≥ 128 mg/L, R) Cefuroxime (≥ 64 mg/L, R) Ceftriaxone (≥ 64 mg/L, R) Cefepime (≥ 32 mg/L, R) Ertapenem (≥ 8 mg/L, R) Imipenem (4 mg/L, R) Meropenem (≥ 16 mg/L, R)	JAWVAI000000000

			Amikacin (32 mg/L, I) Gentamicin (≥ 16 mg/L, R) Ciprofloxacin (> 4 mg/L, R) Fosfomycin (≤ 16 mg/L, S) Colistin (1 mg/L, I) Trimethoprim/Sulfamethoxazole (≥ 320 mg/L, R)	
AMRIR00637	<i>Pseudomonas aeruginosa</i>	<i>bla</i> _{KPC-2}	Piperacillin/Tazobactam (≥ 128 mg/L, R) Ceftazidime (≥ 64 mg/L, R) Cefepime (≥ 32 mg/L, R) Aztreonam (≥ 64 mg/L, R) Imipenem (≥ 16 mg/L, R) Meropenem (≥ 16 mg/L, R) Amikacin (≥ 64 mg/L, R) Ciprofloxacin (≥ 4 mg/L, R) Levofloxacin (≥ 8 mg/L, R) Colistin (1 mg/L, I)	JBLAVS000000000
AMRIR001225	<i>Acinetobacter baumannii</i>	<i>bla</i> _{GES-5}	Piperacillin/Tazobactam (≥ 128 mg/L, R) Ceftazidime (≥ 64 mg/L, R) Cefepime (≥ 32 mg/L, R) Imipenem (≥ 16 mg/L, R) Meropenem (≥ 16 mg/L, R) Amikacin (≥ 64 mg/L, R) Gentamicin (≥ 16 mg/L, R) Ciprofloxacin (> 4 mg/L, R) Colistin (16 mg/L, R) Trimethoprim/Sulfamethoxazole (≥ 320 mg/L, R)	Sanger sequencing
AMRIR00139	<i>Pseudomonas</i>	<i>bla</i> _{GES-9}	Ticarcillin/Clavulanic acid (≥ 128 mg/L, R)	JBBBSP000000000

	<i>aeruginosa</i>		Ceftazidime (≥ 64 mg/L, R) Cefepime (≥ 64 mg/L, R) Aztreonam (≥ 64 mg/L, R) Doripenem (≥ 8 mg/L, R) Imipenem (≥ 16 mg/L, R) Meropenem (≥ 16 mg/L, R) Amikacin (≥ 64 mg/L, R) Gentamicin (≥ 16 mg/L, R) Ciprofloxacin (≥ 4 mg/L, R) Levofloxacin (≥ 8 mg/L, R) Colistin (2 mg/L, I)	
--	-------------------	--	--	--